

# Single-cell RNA-sequencing (scRNA-seq): knowing the in and outs of the data generated

Ken Lau, Assistant Professor of Cell and Developmental Biology

[ken.s.lau@vanderbilt.edu](mailto:ken.s.lau@vanderbilt.edu)

CQS Summer Academy (8/13/2018)

<http://www.mc.vanderbilt.edu/vumcdept/cellbio/laulab/index.html>

Twitter: @KenLauLab



## About me and single cells

- Started lab at Vanderbilt in 2013 with focus on single-cell biology of the gut (IBD and colon cancer)
- Multiplex imaging, CyTOF, scRNA-seq
- Training at Toronto/MIT/Harvard on multivariate analysis, mathematical modeling, and tissue systems
- inDrop in lab since August 2016 (first 1cell customer outside of Boston), > 50 samples ran so far > 500 000 cells sequenced; we have two systems



# Outline

- Introduction to scRNA-seq techniques
- Discussion on scRNA-seq data issues
- Brief Python introduction

[https://github.com/KenLauLab/Discovery\\_Oriented\\_Data\\_Science](https://github.com/KenLauLab/Discovery_Oriented_Data_Science)

## **scRNA-seq techniques**

# Bulk data versus single-cell data






Bulk

VS.




Single-cell

Bulk RNA-seq



1 sample x 25,000 Genes

# Single-cell RNA-seq



# scRNAseq protocols

## Cell Encapsulation techniques

- Droplet-Based
- Well-Based
- Microfluidic capture (Fluidigm C1)

## Lysis and RT

- Coupled – requires balanced mix
- Uncoupled – enables more aggressive lysis

## RNA capture strategies

- Poly dT priming
- Targeting / enrichment

## Indexing strategies


- During capture/RT (typically per **cell** indexes - barcodes)
- After RT (typically per **well** indexes)

## Amplification strategies


- PCR vs IVT





# Typical logic for a scRNA-seq experiment




Single cell suspension\*\*\*



Single-cell encapsulation/  
Library preparation



Sequencing and alignment  
(Bioinformatics I)




# Well-known methods to isolate single cells

## Laser-capture microdissection


(i) Cell selection




(ii) Laser sectioning





(iii) Cell transfer on a membrane



## Fluorescence-activated cell sorting



## Plate-Based scRNAseq





- Isolate RNA, label transcripts using barcoded RT primers (3' seq) or through template switching library prep (enables full length)





Relative to other platforms:

- $\sim 10\mu\text{l}$ / cell,  $<1000$  cells (higher volume, lower throughput)
- Deeper sequencing possible - flexibility
- Number of wells become limiting (doublet rate vs. cost)

# Microfluidic capture scRNAseq (Fluidigm C1)



Capture chips are built to accommodate specific ranges of cell size




Cells

Stanford course

## Bead-based capture: Immobilized sets of indexed primers

- Each bead is coated with primers containing a barcode unique to that bead – (index for each cell)



Barcode



UMI

TTTT

- And a Unique Molecular Identifier (UMI) that uniquely tags each primer – (index for each transcript)

# Droplet-based encapsulation

- Co-encapsulating cells and beads in thousands of 1-5nL droplets
- Beads carry barcoded poly-T primers to capture RNA
- Encapsulation rate follows Poisson distribution
- Excess of “vessels” to minimize doublets



Cells entering/ sec x


Beads/droplet

=

cells captured/ sec

# Beads


- inDrop  
(1cellBio)



- Drop-seq  
(Chemgenes)




- Chromium  
(10x Genomics)





# Microfluidic droplet encapsulation chip

Beads




Cells


Droplet  
formation

Emulsion output


# Seq-Well - microwell sequencing (Shalek lab)




## Drop-Seq




## Seq-Well



# Seq-Well: Principle




## Nanowell Array



**Size Exclusion** →  $\leq 1$  bead per well

**Sealing** → ~~Cross Contamination~~

# Deconvolving the data



Reads with same barcodes collapse into cells

Read with same UMIs collapse into transcript counts

	Cell:	1	2	...	N
Count unique UMIs for each gene in each cell	<i>GENE 1</i>	1	2		14
	<i>GENE 2</i>	4	27		8
	<i>GENE 3</i>	0	0		1
	⋮	⋮	⋮		⋮
	⋮	⋮	⋮		⋮
Create digital expression matrix	<i>GENE M</i>	6	2		0

# **scRNA-seq Data Exploration (and problems)**







## Barcodes


## Genes

0000	0001	0002	0003	0004	0005	0006	0007	0008	0009	0010	0011	0012	0013	0014	0015	0016	0017	0018	0019	0020	0021	0022	0023	0024	0025	0026	0027	0028	0029	0030	0031	0032	0033	0034	0035	0036	0037	0038	0039	0040	0041	0042	0043	0044	0045	0046	0047	0048	0049	0050	0051	0052	0053	0054	0055	0056	0057	0058	0059	0060	0061	0062	0063	0064	0065	0066	0067	0068	0069	0070	0071	0072	0073	0074	0075	0076	0077	0078	0079	0080	0081	0082	0083	0084	0085	0086	0087	0088	0089	0090	0091	0092	0093	0094	0095	0096	0097	0098	0099
0100	0101	0102	0103	0104	0105	0106	0107	0108	0109	0110	0111	0112	0113	0114	0115	0116	0117	0118	0119	0120	0121	0122	0123	0124	0125	0126	0127	0128	0129	0130	0131	0132	0133	0134	0135	0136	0137	0138	0139	0140	0141	0142	0143	0144	0145	0146	0147	0148	0149	0150	0151	0152	0153	0154	0155	0156	0157	0158	0159	0160	0161	0162	0163	0164	0165	0166	0167	0168	0169	0170	0171	0172	0173	0174	0175	0176	0177	0178	0179	0180	0181	0182	0183	0184	0185	0186	0187	0188	0189	0190	0191	0192	0193	0194	0195	0196	0197	0198	0199
0200	0201	0202	0203	0204	0205	0206	0207	0208	0209	0210	0211	0212	0213	0214	0215	0216	0217	0218	0219	0220	0221	0222	0223	0224	0225	0226	0227	0228	0229	0230	0231	0232	0233	0234	0235	0236	0237	0238	0239	0240	0241	0242	0243	0244	0245	0246	0247	0248	0249	0250	0251	0252	0253	0254	0255	0256	0257	0258	0259	0260	0261	0262	0263	0264	0265	0266	0267	0268	0269	0270	0271	0272	0273	0274	0275	0276	0277	0278	0279	0280	0281	0282	0283	0284	0285	0286	0287	0288	0289	0290	0291	0292	0293	0294	0295	0296	0297	0298	0299
0300	0301	0302	0303	0304	0305	0306	0307	0308	0309	0310	0311	0312	0313	0314	0315	0316	0317	0318	0319	0320	0321	0322	0323	0324	0325	0326	0327	0328	0329	0330	0331	0332	0333	0334	0335	0336	0337	0338	0339	0340	0341	0342	0343	0344	0345	0346	0347	0348	0349	0350	0351	0352	0353	0354	0355	0356	0357	0358	0359	0360	0361	0362	0363	0364	0365	0366	0367	0368	0369	0370	0371	0372	0373	0374	0375	0376	0377	0378	0379	0380	0381	0382	0383	0384	0385	0386	0387	0388	0389	0390	0391	0392	0393	0394	0395	0396	0397	0398	0399
0400	0401	0402	0403	0404	0405	0406	0407	0408	0409	0410	0411	0412	0413	0414	0415	0416	0417	0418	0419	0420	0421	0422	0423	0424	0425	0426	0427	0428	0429	0430	0431	0432	0433	0434	0435	0436	0437	0438	0439	0440	0441	0442	0443	0444	0445	0446	0447	0448	0449	0450	0451	0452																																															

# Inflection point method for identifying barcodes with real cells



# Doublet rate determination






# Procedure for the isolation of high viability single-cells from tissues

Keep all reagents on ice and perform all procedures at 4°C, avoid working with overconcentrated tissue/cell solutions


## Tissue

\*Keep cold and minimize ischemic time

\*Mince if needed



Normal colonic epithelium can be isolated through chelation, while other tissues may require direct mechanical processing to achieve 50-500µm fragments. Care should be taken to remove dead cells during washes to maintain viability



The process of isolating tissue fragments should be optimized to the needs of the target cells. This may mean filtering or taking other measures to enrich.

Wash in PBS  
Decant small debris

## Chelation

10ml DPBS (-Mg/Ca)  
3mM EDTA  
0.5mM DTT  
10mM NAC



~30-60min  
Change buffer every 10-15 min


Use forceps to move tissue

Move back to fresh chelation buffer

## Shaking

10ml DPBS (-Mg/Ca) in a new tube

2-3 Shakes/sec  
~6 inches motion




Wash & inspect crypts

## Medimachine

In cold room

Pre-wet 50µm medicon before use  
Wash quickly and thoroughly after




\*GentleMACS is a potentially less damaging approach we have not yet tested for breaking apart tissue and further dissociating cells.

Wash in PBS  
Decant small debris

## Washing


DPBS (-Mg/Ca) spin @~300xg for 2.5min



Good pelleting:  
(fat or pancreas may still float in first wash)



Sticking to tube;  
may need more chelation (EDTA)\*



Floating/clumping  
typically caused by excess mucin or inviable cells\*\*

## Cold Protease Dissociation


2-5ml DPBS (-Mg/Ca) + 2.5mg/ml DNase + 5mg/ml Subtilisin

\*Use 2ml for every 50-100µl of pelleted tissue, don't overconcentrate

\*Add buffer to frozen subtilisin aliquot to thaw immediately before use

\*Pipette with a 1000µl tip every 5-10 min and check for singlets

4-6°C with  
gentle motion



\*Pellet without EDTA prior to re-suspending in protease

\*\*Optimize spin conditions so that single cells are decanted


\*\*Don't leave tissue in pellet too long, split if too dense

\*\*Re-suspend well, but do not over pipette prior to dissociation


# Single-cell isolation for scRNA-seq to minimize dissociation artifacts

- cold protease from *Bacillus licheniformis*, soil bacteria from Himalayan glaciers
- enables tissue preparation on ice (at 4 degrees)

37 degree DNase/Collagenase




Cold Protease (on ice)




- caveat – the efficacy for retrieving all cell types from all tissues unknown

# Stressed/dying cells


PCA of mitochondrial gene expression



# Contaminant from free floating RNA and leaky cell corpses




Generous threshold



Strict threshold


# Amplification bias in scRNA-seq




Non-linear amplification

- Highly expressed transcripts are inappropriately represented and replicated
- Sampling of RNA in a cell 1-10%
- Zero inflated data
- Count data – negative binomial distribution

# In vitro transcription results in linear amplification of RNA




# Negative binomial versus Gaussian distributions



## Gaussian

- Continuous variable
- Symmetric



## Negative Binomial

- Discrete variable (counts)
- Asymmetric at small means



# Zero inflation in scRNA-seq

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25
1	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
2	0	0	0	0	0	0	0	1	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
3	0	0	0	0	0	0	0	0	0	0	2	0	0	0	0	0	0	0	0	0	0	0	0	0	0
4	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	1	0	0	0	0
5	0	0	0	0	0	0	0	0	0	0	2	0	0	0	0	0	0	0	0	0	0	0	0	0	0
6	0	0	0	0	0	0	0	1	0	0	1	0	0	0	0	0	1	0	0	0	0	0	0	0	0
7	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	1	0	0	0	0
8	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
9	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
10	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
11	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
12	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	1	0	0	0	0
13	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	1
14	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
15	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	1
16	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	1	0	0	0	0	0	0	0	0
17	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
18	0	0	0	0	0	0	0	2	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
19	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
20	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
21	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
22	0	0	0	0	0	0	0	1	0	0	1	0	0	0	0	0	0	0	0	0	0	0	0	0	0
23	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
24	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
25	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
26	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	2	0	0	0	0
27	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
28	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
29	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
30	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
31	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
32	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	1
33	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	1	0	0	0	0	0	0	0	0
34	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
35	0	0	0	0	0	0	0	0	0	0	2	0	0	0	0	0	0	0	0	0	0	0	0	0	0
36	0	0	0	0	0	0	0	0	0	0	1	0	0	0	0	0	0	0	0	0	0	0	0	0	0
37	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0

- Only 5% of this table is not 0!



# What does this mean?

- Statistical tests to determine differentially expressed genes with a zero-inflated negative binomial (ZINB)
- False negatives (genes that are supposed to be expressed by appear as 0)
- Many genes (columns) have low counts due to shallow sampling of transcripts – this means the data are noisy
- Unreliable variables that need to be processed/filtered out prior to downstream analysis

# **Feature Selection**

# Variance selection


- Easiest – select genes that are the most variable

- Variance = genes that are most different across all cells

$$\sigma^2 = \sum \frac{(X - \mu)^2}{N}$$

- Rank genes by top 500 most variable, for example, discard the rest


# Highly variable axes = Principal Components



# PCA (Principal Component Analysis)

- Principle of PCA is to maximize the Variance of  $X$  with the least amount of principal components (latent variables)
- What is variance? Spread of the data, information content, change etc.
- Variance is the covariance of a dataset with itself, i.e.  $\text{Var}(X) = \text{Cov}(X, X) \rightarrow \text{Maximize}$
- What are principal components? Linear combinations of original variables – linear transformation


# Vectors and projections



Basis Set

$(1,0)$   $\longrightarrow$   
 $(0,1)$   $\uparrow$


$\longrightarrow$   
Projection




Basis Set

$(7,4)$   
 $(-3.5,7)$

# PCA as a dimension reduction tool



# Selection of principal components



**PRINCIPAL COMPONENTS CAN BE ORDERED BY EIGENVALUES (VARIANCE CAPTURED)**



## Simple summary of this simple feature selection procedure

- Keep most variable genes (over all cells) for downstream analysis, discard rest as noise
- PCA identifies super axes (Principal Components) that are combinations of the original variables that captures the most variance in the data (by eigenvalues)
- Orthogonal – no duplicate or redundant axes – so will only have a few of them